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Operation Instructions

General Sonic® GS2 – GS4

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www.qteck.de

Ultrasonic cleaners for aqueous cleaning liquids

Read this instructions before you start operating your ultrasonic cleaner and keep this booklet at hand for all persons that might use this device.

For more information about your unit and other qteck products please visit us at:

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Many thanks that you have decided for a **General Sonic** ultrasonic-cleaning device. These operating instructions should transmit to you important and helpful tips and safety regulations in working with your new ultrasonic-cleaner. Into relation on the choice of cleaning concentrate, cleaning temperature and duration of insonification can be entered here on grounds of the huge number of the possibilities which arise with the help of the big operational areas not on all details.

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Table of contents:

indication	page
1. General safety and operating information	9
2. Commisioning	9
3. Ultrasonic cleaning process	10
3.1 Principle of ultrasonic cleaning.....	10
3.2 Helpful tips	10
3.3 Kinds of cleaning	10
3.4 Operating elements.....	11
3.5 Cleaning procedure	11
3.6 After the cleaning process	12
3.7 Preparation of medical instruments	12
4. How to clean the device	12
5. Funktional disorder.....	12
6. Technical specification.....	13
7. Warranty	13
8. Declaration of confirmity.....	14

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1. General safety and operating information

- Devices, accessories and cleaning agent may be used only according to the operating instructions or the manufacturer's data.
- Ultrasonic-cleaning devices are not made for long-term irradiation (permanently pursued) and need short quiescence between the cleaning processes
- The device may be pursued only with the mains voltage given on the type sign (device back) in an earthed protective contact outlet
- Devices with damages may not be connected to the power supply system.
- Devices and cleaning agents are to be kept away from children and not authorised persons
- For material damages or personal damages that are caused by improper use or nonobservance of the security terms, no liability is taken over
- By use of cleaning agents the security and application instructions of the manufacturer must be exactly considered
- High-handed doing alterations or changing the devices is not permitted, in addition, the CE-conformity expires.
- Use exclusively admitted accessories only.
- In commercial facilities the accident prevention regulations of the federation of the commercial professional association are to be followed for electric arrangements and company means.
- In schools, training centres, amateur workshops and self-help workshops the incentive is to be supervised by trained staff responsibly.
- An auditory protection is to be carried (worn) during the company (operation) for the protection of the health on constant (continuing) activity in the vicinity of less than 3 m to the device (appliance).
- Do not pursue device without or with liquid not enough.
- While insonification of small quantities of ignitable liquids in insert vessels the guidelines are to be followed for laboratories BGR120 and the explosion protection guidelines (EX-RL / BGR104) in the each case valid version.

2. Commissioning

- Put the device on a steady, horizontal and flat surface.
- Fill in the tank to 3/4 with water. **Note:** The liquid temperature may not be higher than 50°C, never pursue devices without liquid.
- Add then the cleaning agent (e.g., CQ55, CQ11, CQ33) according to the information of the manufacturer in the right pouring proportion. **Note:** use only admitted cleaning agents. No ignitable liquids (e.g., petrol, solvent) and no liquids which contain chloride ions or split off (some disinfection, washing-up liquid or budgetary cleaner) use. Aggressive cleaning liquids as for example acids or salt solutions are allowed only in closed insert vessels.
- Put the power plug in an enough secure outlet.

Tip:

Only after approx. 1-2 hours of operation the transducers (sonic generators) works with full power. Cleaning liquids and contact liquids contain relaxed gases (e.g., oxygen) that the cleansing effect negatively influences. Degas freshly filled liquid before use by insonification approx. for 10 min. Last to this process the cavitation noise becomes quieter, the loud degassing noises are cancelled and the device works noticeably more quietly. Nevertheless, this signifies nobody decrease to the ultrasonic achievement separate points the end of the degassing process.

3. Ultrasonic cleaning process

3.1 Principle of ultrasonic cleaning

The transducers fixed under the liquid tank convert the electric energy into mechanical oscillations. Then the cleaning liquid is brought with 35 or 45 kHz (according to device) to swing. Besides, originates a cavitation with which small (little) vacuum vesicles implode and are thereby removed from the deepest pores pollutions. One calls (mentions) this process also „ electronic brushes “.

3.2 Helpful tips

- Use only special ultrasonic cleaners (e.g., CQ11, CQ33 or CQ55) as an addition. Tap water without any addition cleans insufficiently. "Tip degasing follow ".
- Warm liquids raise the cleansing effect and shorten the cleaning duration.
- Heat the liquid before the insertion of the cleaning property on the desired temperature. Pay attention, on this occasion, absolutely to the filling state of the liquid, with too low filling state damages can appear in the heating. To save energy and time, you can cover the liquid tank with the lid.
- At too high temperatures the ultrasonic - cavitation decreases again, the empirical best cleaning temperature lies with approx. 50 - 60°C.
- The cavitation warms up the liquid, in addition. By long duration of ultrasonic irradiation and/or by covering the tank, the temperature of the liquid can also increase above the set value of the thermostat.
- Disinfection liquids may not be warmed up, in addition.
- Make sure that the liquide level is constand during the whole cleaning process.
- Strongly dirty parts preclean mechanically and insert with the stronger dirty side down. Do not pile.
- In hollow cavities no air bubbles may be

3.3 Kinds of cleaning

- 1. Direct cleaning in the swinging tub / tank**
- 2. Indirect cleaning in closed application vessels**

3.3.1 Direct cleaning


The ordinarily applied kind of the cleaning occurs directly in the liquid located (contained) in the tank. Moreover the cleaning object is brought in an basket or a device holder in the prepared liquid. Basically is to be paid attention to the fact that the cleaning object is covered completely with liquid. Not overcast places are not cleaned. No accessories (basket) or cleaning goods may touch the tank ground. Put bits and pieces in an additional sieve basket.

3.3.2 Indirect cleaning

An indirect cleaning in closed application vessels is carried out with:

- Application from chemically aggressive (e.g., acids) or ignitable (e.g., petrol) cleaning liquids.
- Distance of chemically aggressive soiling.
- Distance of abrasive soiling (e.g., polishing pastes, quartz, sand).
- concurrent use of different cleansing liquids

3.4 Operating elements

- 1 Power switch.....switch on and off the Unit
- 2 Timer switchto adjust time of procession / start ultrasonic
- 3 Pilot lamp.....shows heating works
- 4 Temperature switchto adjust heating temperature

3.5 Cleaning procedure

- Switch on the device at the power switch.
- Turn the temperature selector switch on the desired cleaning temperature. (the controlling light (3) begins with some delay), the heating works.
- After the cleaning temperature is reached lay the cleaning objects in the basket and hang up this carefully in the tank. Attention: There is the danger of scald by splash, in the liquid do not touch. Cleaning objects turn around perhaps with tools to avoid air bubbles inside. With the bringing in of the basket pay attention to the fact that the tank is not too full, or overruns. In that case, liquid must be taken out accordingly.
- Turn up the timer for the maximum time and then back for the desired insonification time.
Tip: Choose the insonification time so briefly as possible. According to soiling, size and used cleaning liquid the cleaning duration amounts experience to approx. 5 - 15 min.
- With parts (sections) sensitive to warmth follow temperature (distortion danger).
- At the end of the insonification time take out the cleaning objects with the basket (care combustion danger). Let the cleaning objects cool down and rinse it under fluent water.
- If necessary repeat cleaning process
- With the indirect cleaning the cleaning objekt is laid(put) in an e.g. glass application (as accessories available) and is covered with cleaning liquid and put(arranged) in the basket. (or with the matching hold ring in the hole lid (as accessories available, for 2 glass applications, for the use of different cleansing liquids) hung up. **Attention:** overcrowding danger.
- The application vessels must dip min 3 cm deep into the contact liquid.
- You can use as a tank filling either the available (in existence) cleaning liquid or special contact liquid for ultrasonic – transference (transmission).

3.6 After the cleaning process

After every cleaning process the tank should be cleaned. On the ground of the tank mature dirty leftovers reduce the ultrasonic achievement and can lead to damages on the tank ground. Cleaning liquids can be used according to soiling several times, but they cannot be refreshed by recent addition of cleaner again. Spent cleaning liquids have a lower cleaning effect and have to exchange.

3.7 Preparation of medical instruments

- If use the ultrasonic - cleaner in connection with disinfection preparations (medicine product class IIa) in addition to support the disinfection effect or to accelerate as accessories to a medicine product and is to be arranged even as a medicine product class I.
- In case of cleaning or disinfection of medical instruments, the hygienic security is to be guaranteed after the utilisation. The cleaning and disinfection has to be carried out by the operator, according to the hygiene guideline (plan) under using a e.g. VAH/DGHM certified sureface disinfectant.

4. How to clean the device

- Before cleaning every device disconnect from power supply system.
- By devices without draining cock pour out the tank contents over the corner (backside left), besides, splashes avoid to the device ground or the cable connection.
- Devices do not wash up or dip into water.
- Rinse the tank thoroughly and rub dry, no steel wool, scratches, scrapers use or similar.
- Edges and remains in the tank are to be removed with a customary high-grade steel cleaner without abrasive addition. **Tip:** Remaining metal parts on the stainless steel surface as well as rust particle from the water pipe system or the cleaning property can penetrate the passive protective coating of the stainless steel, "activate" the stainless steel and it starts to rust. This foreign rust causes pitted areas and the tank becomes leaking.
- Wipe case outside to the cleaning only humid, use when required a customary high-grade steel cleaner without abrasive addition.

5. Functional disorder

problem	possible cause	remedy
weak uneven insonification, noise is too loudly or uneven	liquid is not degased	10 min. insonification
	too mutch cleaning propertys	put some out
	uneven noise	change filling state
Liquid does not become warm	filling state to low	fill up liquide
	heating damaged	customer service
pilot lamp does not shine	lamp or heating damaged	customer service
liquid to hot	warm by kavitation	put down temperature
	tank covered	remove lid
Device cannot be switched on	electricity supply discontinuous	check, produce elektricity supply
no insonification	electronic problem	customer service

6. Technical specification

Device type	GS-2	GS-3	GS-4
Tank volume ml	2000	3000	4000
Power input Watt	250	450	500
Maximum power	360	600	653
Protection class	IP 20	IP 20	IP 20
Maximum Ultrasonic power HF Watt	153 140	384 210	360 280
Ultrasonic transducers	2	3	4
Frequency kHz	35	35	35
Heating adjustable -°C	20 – 80	20 - 80	20 - 80
Maximum Heating Watt	230 150	300 300	550 300
Timer min.	0 – 30	0 – 30	0 – 30
Case	stainless steel	stainless steel	stainless steel
Case measures mm	260 x 150 x 210	260 x 160 x 228	323 x 177 x 228
Tank measures mm	213 x 99 x 53	237 x 137 x 99	298 x 154 x 99

6.1 Information for use as a medical product


Denotation	ultrasonic cleaning device
UMDNS - nomenclature (ECRI)	14-263
Purpose determination	ultrasonic-disinfection and cleaning of medical and dental instruments
Classification (Medical device directive 93/42/EEC, appendix IX)	class I, active, non invasive, not implantable medical product
Device specification	see type plate on device backside

7. Warranty

The warranty duration for material defect liability on our devices is 24 months for private users and 12 months for commercial users. The warranty duration begins with the purchase of the device. Defects are to be indicated immediately by writing statement and the device must be put except work. In the warranty case the manufacturer reserves himself, in any case, the right on finishing touches. An occurred exchange, whether on the whole or partially extended not the warranty duration. For the damages on the device or cleaning objects which are caused by improper use, use of inexpedient cleaners or chemicals no warranty exists.

Many damages are mostly to be led back on operating mistakes. Please, note:

- **Damages in electronics / ultrasonic generators = too long irradiation times / the longest allowed time in one piece is 30 min..**
- **Distortion of the cleaning objects = too high temperature**
- **Heating defect = filling state too low, purposed with not enough liquide**
- **corrosion (rust) in the tank = bad, insufficient cleaning of the device, wrong cleaners, aggressive liquid directly applied.**
- **Damages by liquid in the electronics = tank leaking, humidity from the outside, penetrated in by device cleaning.**

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