Product Description

3M Filtek Easy Match Universal Restorative material, is a visible-light activated, radiopaque restorative composite designed for use in anterior and posterior restorations. The three available shades (Bright, Natural, Warm) have a body-like opacity. The fillers are a combination of a non-agglomerated/non-aggregated 20 nm silica filler, a non-agglomerated/non-aggregated 4 to 11 nm zirconia filler and an aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles). 3M Filtek Easy Match Universal Restorative has an average cluster particle size of 0.6 to 10 microns. The inorganic filler loading is 78.5% by wt. (63.3% by volume). 3M Filtek Easy Match universal contains bis-GMA, UDMA, TEGDMA, PEGDMA and bis-EMA resins. A dental adhesive, such as those manufactured by 3M Company, is used to permanently bond the restoration to the tooth structure. It is packaged in traditional syringes and single-dose capsules.

General Information

3M Filtek Easy Match Universal Restorative complies with ISO 4049: Polymer-Based Dental Restorative Material classified as a Type 1 and Class 2 material.

All shades are radiopaque, with a value of 2.4 mm of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1 mm of material having a radio-opacity equivalent to 1 mm of aluminum has a radio-opacity equivalent to that of dentin and 2 mm of aluminum is equivalent to enamel.

Intended Use

A dental composite resin for anterior and posterior restorations.

Indications for Use

3M Filtek Easy Match universal restorative is indicated for use in:

- Direct anterior and posterior restorations (including occlusal surfaces)
- Core Build-ups

Contraindications

None

Intended Users

Educated dental professionals, i.e., general dentists, dental assistants/hygienists, who have theoretical and practical knowledge on usage of dental products.

Intended Patient Population

Patient requiring dental restoration as recommended by a dentist unless the patient's condition, such as a known allergy to the device, limits the use.

Clinical benefit

Restoration of oral esthetics and function.

Risks and Undesirable effects

- Allergic reaction (i.e., redness, swelling, itching, burning sensation near the site of application)
- Post-operative sensitivity/pain

Precautions and Warnings For Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

For Dental Personnel

Capsules may be warmed (Do not warm syringes).

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If an allergic reaction occurs, seek medical attention as needed.

Instructions for Use

Preparation

1. Prophylaxis: Teeth should be cleaned with pumice and water to remove surface stains.

2. Shade Selection: Before isolating the tooth, select the appropriate shade of restorative material. In most restorations, the shade "Natural" will match. For brighter or warmer smiles, the "Bright" or "Warm" shade can be selected.

2.1 Initial product use / tooth shade is difficult to identify:

- a) Button matching: Observe the patient's tooth to be restored. Select "Natural" shade, place a small button of composite on the tooth to be restored (without etching or bonding). Manipulate the material to approximate the thickness and site of the restoration. Cure. Evaluate the shade match under different lighting sources.
- Good shade match, proceed with selected "Natural" shade
 - If the composite is too light compared to the tooth, select "Warm" shade
 - If the composite is too dark compared to the tooth, select shade "Bright" shade

Remove the restorative material from the unetched tooth with an explorer.

b) Conventional approach: Determine the shade of the tooth to be restored by using the VITA classical shade guide. Select the corresponding shade from the table below.

Bright	Natural	Warm
Bleach, A1, B1, B2, C1, D2	A2, A3, C2, D3, D4	A3.5, A4, B3, C3, C4

2.2 Intuitive product use: Observe the patient's tooth to be restored. Intuitively select the shade depending on the chroma and value. In most restorations, the shade "Natural" will match. For brighter or warmer smiles, the "Bright" or "Warm" shade can be selected.

Note: For complex multi-shade restorations, consider a multi-opacity, multi-shade composite system.

3. Isolation: A rubber dam is the preferred method of isolation. Cotton rolls plus an evacuator can also be used.

Direct Restorations

1. Cavity Preparation:

1.1 Anterior restorations: Use conventional cavity preparations for all Class III, IV and V restorations.

1.2 Posterior restorations: Prepare the cavity. Line and point angles should be rounded. No residual amalgam or other base material should

be left in the internal form of the preparation that would interfere with light transmission and therefore, the hardening of the restorative material.

2. Pulp Protection: If a pulp exposure has occurred and if the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of a Liner/ Base as appropriate, such as those manufactured by 3M Company.

3. Placement of Matrix: Place the matrix system of choice by following the manufacturer's instructions for use.

4. Adhesive System: Follow the manufacturer's instructions regarding etching, priming, adhesive application, and curing, for example 3M™ adhesives.

5. Dispensing the Composite: Follow the directions corresponding to the dispensing system chosen.

5.1 Syringe: Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of the restorative when dispensing is completed, turn the handle counterclockwise a half turn to stop paste flow. Immediately replace the syringe cap. If not used immediately, the dispensed material should be protected from light.

5.2 Single-Dose Capsule: Insert capsule into 3M™ Restorative Dispenser, manufactured for 3M Company. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

6. Placement:

6.1 Place and light cure restorative in increments, as indicated in Section 7.

6.2 Avoid intense light in the working field.

6.3 Posterior placement hints:

a) To aid in adaptation, the first 1 mm layer may be placed and adapted to the proximal box.

b) A condensing instrument (or similar device) can be used to adapt the material to all internal cavity aspects.

7. Curing: This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 550 mW/cm² in the 400-500 nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M™ curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

Increment depth	Cure Time	
	All halogen or LED lights LED lights (with output 550-1000 mW/cm ²)	Elipar™ S10 (3M™ LED lights with output 1000-2000 mW/cm ²)
2.0 mm	20 sec.	10 sec.

8. Contouring: Contour restoration surfaces with fine finishing diamonds, burs, or stones. Contour proximal surfaces with Sof-Lex™ Finishing Strips, manufactured for 3M Company.

9. Adjust Occlusion: Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

10. Finish and Polishing: Polish with the Sof-Lex™ Finishing and Polishing System.

Storage Conditions

1. This product is designed to be used at room temperature. If desired, the product may be warmed in a commercial warmer prior to use (no higher than 70C/158F, no longer than 1 hour); for capsules only.

2. The product is best stored at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 36 months. Ambient temperatures routinely higher than 27C/80F may reduce shelf life. See outer package for expiration date.

3. Do not expose restorative materials to intense light.

4. Do not store materials in proximity to eugenol containing products.

Handling after Use

Capsules

The capsule is single-use only and should be discarded after use.

Syringes

Handling

To prevent contamination, the syringe should be handled outside the treatment area with new, unused disposable gloves. The syringe cannot be reprocessed. Dispose the syringe if contaminated, e.g., with blood, saliva, spray of body fluids, or contaminated hands. The syringe will not be damaged by incidental contact with water, soap, or standard cleaning and disinfecting wipes or spray. Repeated liquid contact and vigorous wiping may damage the label. Do not put the syringe into a disinfection bath or a washer-disinfector (thermal disinfector).

Inspection, Maintenance, Testing

Examine the syringe for damage, discoloration, and contamination before every use. Do not use damaged and contaminated syringes under any circumstances.

Ingredients

Silane Treated Ceramic, Bisphenol A Diglycidyl Ether Dimethacrylate, Bisphenol A Polyethylene Glycol Diether Dimethacrylate, Diurethane Dimethacrylate, Silane Treated Silica, Polyethylene Glycol Dimethacrylate, Silane Treated Zirconia, Triethylene Glycol Dimethacrylate


SDS and Disposal

- See the Safety Data Sheet (available at www.3M.com) for disposal and ingredient information.
- See SDS for disposal information or to be disposed in accordance with the local or regional regulations.
- Please pay special attention to the disposal of contaminated waste to avoid health risks due to improper handling.

Customer Information

- Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.
- No person is authorized to provide any information which deviates from the information provided in this instruction sheet.
- A Summary of the Safety and Clinical Performance can be found at <https://ec.europa.eu/tools/eudamed>, Basic UDI-DI: 06082238401020000000047BK

3M and Filtek are trademarks of 3M. All other trademarks herein are the property of their respective owners.

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Non-corrugated fiberboard		Indicates cardboard material is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Temperature Limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Refer to instruction manual/booklet		To signify that the instruction manual/booklet must be read. ISO 7010-M002
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU

For more information see, HCBGregulatory.3M.com

(*) This symbol is blue on the packaging

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3M, Filtek and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. © 3M 2023. All rights reserved.
3M Deutschland GmbH
3M (Schweiz) GmbH